Tizian Blank translucent

CE 0297

Schütz Dental GmbH · Dieselstr. 5-6 · 61191 Rosbach/Germany
Telefon: +49 (0) 6003 814-0 · Telefax: +49 (0) 6003 814-906
info@schuetz-dental.de · www.schuetz-dental.de

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Instructions

Description / indications

Tizian Blank translucent is made of a highly translucent zirconium dioxide material and is indicated for the production of full zirconium bridges. Therefore, Tizian Blank translucent is eminently suitable for coloring even without a ceramic veneer. Tizian Blank translucent is particularly applicable in case of space constraints, for implant-borne or for gingival-borne prostheses. Tizian Blank translucent can be designed in a fully anatomical fashion as well as in the classic technique.

Indications

The blanks are suited to the production of:

- crown copings for anterior and posterior restorations (fully anatomical or reduced)
- fully anatomical anterior or posterior bridges up to 14 units (anterior: max. three connecting links; posterior: max. two connecting links)
- reduced anterior or posterior bridges up to 5 units (anterior: max. two connecting links; posterior: max. two connecting links)
- abutments, inlays
- primary conus crowns and telescoping crowns

Safety note

Tizian Blank translucent is supplied in presintered condition. The blanks show only a limited solidity and exhibit a residual porosity. Therefore, it is necessary to handle them carefully. Please check immediately after receiving the consignment whether packaging and product are intact.

Processing of Tizian Blank translucent

1. Preparation guidelines

Chamfer or shoulder preparations are appropriate for zirconium dioxide. Unsuitable are e.g. channel preparations or tangential preparations. Those preparation types could lead to failure. The following preparation angles are recommended: horizontal -4° to 6° ; vertical -3° to 5° . All connecting areas must be designed in a rounded and smooth fashion. Avoid sharp angles. Provide an optimal insertion vector between the posts. Remove at least 1.5 mm of substance from the oc-

clusal surfaces of posterior teeth and 2 mm from the incisal surfaces of anterior teeth.

2. Milling

Only use machines, systems and tools to process Tizian Blank translucent which are suited to milling presintered zirconium dioxide materials. The blanks must not come into contact with the coolant. Please adhere to the instruction manual for your milling machine.

3. Dimensioning of the framework

After sintering, the wall thickness must be at least 0.7 mm. Finally, after dense-sintering, the crown walls must not be thinner than 0.4 mm. The margins should be at least 0.2 mm thick. After dense-sintering, the connector dimensions in the anterior area should not fall short of 9 mm² and those in the posterior area must be at least 12 mm². The connector geometry is of greatest significance for the crack resistance of zirconium dioxide frameworks. Therefore, choose the greatest connector height possible. Do not place more than two bridge units between the posts in the posterior area. Place only one bridge unit on a lateral free-end. Rule of thumb: The more bridge units between the posts, the more stable must the framework be designed. Design the framework in anatomic fashion to avoid chipping. This way, the veneering ceramic can be applied in regular layers with a maximum of 1–2 mm thickness.

4. Coloring

All conventional zirconium dioxide coloring liquids can be used. Please refer to the manufacturer's instructions.

Obtain best possible coloring results by individually painting the restoration with a brush. Dipping the restoration into the stain is not recommended. The painting technique is similar to the technique used before glaze firing. Use a brush made only of plastic or wood, without any metal parts, as metal parts could cause discoloration. Please observe the instructions of the stain manufacturer.

5. Sintering in a high-temperature furnace

Remove all milling dust from the frameworks before dense-sintering. By sintering, the restoration which was produced about 25% bigger than necessary, will be shrunk to its final size. During sintering, the restoration will obtain its final characteristics. It is absolutely necessary to use a suitable furnace which is approved for presintered zirconium dioxide objects.

Sintering program:

Start: room temperature

Heating rate: 3°C/min. for bridges

4°C/min. for crowns

Final temperature: 1530°C

Holding time at final temperature: 120 min.

Cooling rate: 7°C/min.

Open the furnace at: 65°C

Total time: 14–16 hrs.

Ensure that wide-span constructions are well supported, thereby avoiding sintering distortions. Zirconium dioxide is a bad heat conductor which means that the zirconium dioxide of massive bridge units is significantly less apt to be damaged by heat than that of thin coping walls.

6. Finishing the restorations

Use a handpiece for finishing an unsintered object. Suitable are zirconium dioxide stones, sintering diamonds, rubber polishers or carbide cutters.

If at all possible, dense-sintered objects should not be finished mechanically. If mechanical finishing cannot be avoided, please only use immaculate water-cooled diamond tools to avoid damaging the structure of the object. Avoid extreme local heat (e.g. when using steam jets, sand-blasters or when polishing to high lustre).

Do not apply excessive pressure during finishing and take special care when working with very graceful constructions. Ensure that the wall thickness does fall short of 0.4 mm. Preferably, do not grind the interdental connectors. In principle, do not grind the connectors from basal direction or apply chamfers. Avoid sharp edges and the use of separating discs. Sand-blasting is not recommended. If

sand-blasting is absolutely necessary, however, it should be done with 120 µm aluminium oxide at less than 2 bar. Rinse the object with water afterwards.

If the sintered objects must be ground extensively, carry out a special correction brand to relieve any tensions in the object: Heating rate 8°C/min., final temperature 1055°C, holding time 5 min., natural cooling.

7. Veneer (if desired)

Only use veneering ceramics suited to veneering zirconium dioxide, e.g. Tizian ceramic. Please adhere to the manufacturer's instructions.

To obtain an optimal bond, we suggest a wash firing with dentin mass (100°C higher than the normal firing temperature) in addition to the bonding agent, e.g. Tizian Liner, recommended by the manufacturer.

Zirconium dioxide is a bad heat conductor which means that the zirconium dioxide of massive bridge units is significantly less apt to be damaged by heat than that of thin coping walls. This fact leads to special firing procedure. If there are no massive bridge units in the bridge to be veneered, we recommend to use the normal heating rate of 55°C/min. Should the volume of the bridge units increase, the compensation of the different heat conduction between crowns and bridge units requires increasing the holding time at final temperature to 2 min. It is always a good solution to lower the temperature by about 10°C to keep the ceramic from fusing too strongly. In case of ceramics layered on zirconium dioxide restorations, a 2 min. holding time at final temperature under vacuum should generally be observed.

Attention: The more massive the framework, the lower the necessary heating rate. Likewise, the objects should be cooled down slowly to avoid cracks. a holding time of 2–3 min. must be observed. Absolutely avoid thermal shocks during the firing procedure (especially when firing massive frameworks). Do not remove frameworks from the furnace before they have been cooled down to at

least 200°C. Never place warm objects onto a cold surface (e.g. table top) as this could cause cracks.

8. Fixation

The restoration can be fixed with phosphate or glass ionomer cement, or by adhesion. Only use materials compatible with zirconium dioxide to integrate the restoration. Observe the manufacturer's instructions.

Please note!

When milling partially sintered or dense-sintered zirconium dioxide, respirable dust can be discharged which may damage the lungs or lead to irritations of the eyes and skin. Therefore, zirconium dioxide must only be processed while using a proper and efficient extractor system. Technicians are recommend to wear protective masks and goggles.

Storage

Store in a dry place at temperatures between 10°C and 50°C (50°F and 122°F). Do not subject the blanks to thrusts or vibrations. Generally avoid any contamination.

Disposal

The material can be disposed of with normal household trash.

Guarantee

Our processing instructions, regardless of whether they are issued verbally, in writing or in the form of a practical demonstration, are based upon our own experience and can therefore be regarded solely as guidelines.

Side effects

Undesirable side effects of this medical device are to be expected only in extremely rare cases if the product is handled and used correctly. However, immune reactions (e. g. allergies) or local discomfort cannot be excluded completely. If you become aware of undesirable side effects, please notify us, even in case of doubt.

Contra-indications/interactions

If a patient has known allergies against or hypersensitivities towards a component of this product, we recommend not to use it or to do so only under strict medical supervision. In such cases, we will supply the composition of our medical device upon request. The dentist should consider known interactions and crossreactions of the product with other materials already in the patient's mouth before using the product.

Note

Please ensure that all above information is readily available to the attending dentist in case this product is used. Please also note the material safety data sheets available.

Dental ceramic in accordance with DIN EN ISO 6872, type II, class 6